

# Luminescence Ratio Analyser toolset for Fiji

Yan CHASTAGNIER

November 14, 2016

## Contents

|          |                                |           |
|----------|--------------------------------|-----------|
| <b>1</b> | <b>Setup</b>                   | <b>2</b>  |
| <b>2</b> | <b>Tools description</b>       | <b>2</b>  |
| <b>3</b> | <b>Image naming convention</b> | <b>3</b>  |
| <b>4</b> | <b>Tools use</b>               | <b>4</b>  |
| 4.1      | Clean Tool . . . . .           | 4         |
| 4.2      | Crop Tool . . . . .            | 5         |
| 4.3      | Divide Tool . . . . .          | 6         |
| 4.4      | Analyse Tool . . . . .         | 8         |
| 4.5      | Param Tool . . . . .           | 11        |
| 4.6      | Min Max Tool . . . . .         | 12        |
| <b>5</b> | <b>Output files structure</b>  | <b>12</b> |
| <b>6</b> | <b>Changing default values</b> | <b>13</b> |

# 1 Setup

Put .ijm file into the folder Fiji.app/macros/toolsets/ (on mac os, right click Fiji.app > Show content folder).




Figure 1 – Location to put .ijm file.

Click on » symbol at the extrem right of the toolbar then select the toolset to display the tools.




Figure 2 – Display the tools.

## 2 Tools description

### 2.1 Clean tool

Apply a median filter of radius 1.

Remove the background signal by subtracting the median value of a region that contains only noise.

Align the stacks to correct the small xy drift due to the multipositioning of the microscope by computing the translations coordinates for a stack of images and applying it to both (optional).

## 2.2 Crop tool

Extract sub regions to analyse them separately with different thresholds.

## 2.3 Divide tool

Set to 0 all points below a computed dynamic threshold on the clean donor stack of images. The threshold can be obtained using different methods. Automatic Threshold methods which compute the threshold level based on the whole image (or stack). Automatic Local Threshold methods compute threshold for each pixel based on the surrounding pixels in a given radius. Manual methods: selecting an area and calculating its median or mean \* coefficient.

Divide pixel by pixel the acceptor images by the donor images, to obtain ratiometric images. Set the image to 16-colors and adjust the range of values to visualize it in pseudo-colors.

## 2.4 Analyse tool

Select regions of interest and measure the mean and standard deviation on each slice. Plot data for each regions.

## 2.5 Close tool

Close ROI Manager, Results window, Log window and all images.

## 2.6 Param tool

Change behavior of other tools and set default parameters.

## 2.7 Min Max tool

Allows to change the visualisation range of all "Ratio" images in a folder and its subfolders, in order to have an homogeneous distribution across images.

# 3 Image naming convention

It is suggested to name raw image files xxxDonorName.tif and xxxAcceptorName.tif, where xxx can be anything as long as it's the same for donor and acceptor images. DonorName and AcceptorName have to be entered in the parameters into the fields "Donor" and "Acceptor" respectively. This way, when choosing the donor, acceptor will be automatically found and both will be opened (this also works if DonorName and AcceptorName are in the middle of the image names). This is a necessary step to do batch processing.

It is possible to use any name, but both donor and acceptor will have to be selected, and the subfolder will be the donor name without .tif extension, if "DonorName" is not found in the file name. The DonorName and AcceptorName (which are in the fields "Donor" and "Acceptor" in the parameters) will be used in the names of the clean images.

## 4 Tools use

### 4.1 Clean Tool

1. Select raw donor image (and raw acceptor image if necessary).




Figure 3 – Window to select raw donor image.

2. Select background area then click "OK".

A 32\*32 pixels area minimizing mean value over the stack is automatically selected. If this value seems always good, consider ticking *Batch clean?* in the parameters to speed up the process.




Figure 4 – Background area selection.

3. Choose if you want to align images in the stack or not.




Figure 5 – Choose to align images or not.

4. Images are processed then saved in a subfolder.

Note: If *Batch clean?* is ticked in the parameters, replace all steps by "Select the directory containing raw images".

## 4.2 Crop Tool

1. If images are not already open, select donor image (and acceptor image if necessary).
2. Select a rectangular area to crop, then click "OK".



Figure 6 – Crop area selection.

3. Enter a label for the crop then click "OK" (this step is skipped if *Give a label to CROPS?* is unticked in parameters).




Figure 7 – Give a label to crop.

4. Select nothing then click "OK" to stop making crops.
5. Images are saved then closed.

Note: If *Batch crop?* is ticked in the parameters, replace step 1 by "Select the directory containing raw images", then apply other steps for each pair of images.

### 4.3 Divide Tool

1. If images are not already open, select donor image (and acceptor image if necessary).
2. If threshold method is "Median" or "Mean \* coefficient", select area then click "OK".
3. The threshold is computed and applied for each slice and displayed in a new window. Check if the result is good then click "OK".




Figure 8 – Check threshold.

4. Adjust threshold parameters if needed. Keep *Update Threshold Settings?* ticked to start thresholding process again with new parameters. Untick it to keep settings and continue.




Figure 9 – Threshold settings.

5. Ratiometric image is obtained dividing acceptor by donor. It is displayed in 16 colors. Check the color range (click Live on slice histogram to see distribution along the stack, or use the stack histogram), then click "OK".




Figure 10 – Check color range.

6. Adjust color range parameters if needed. Keep *Update Min and Max?* ticked to update the range. Untick it to keep settings and continue.




Figure 11 – Color range settings.

7. Ratiometric image is saved then closed.

## Notes

- In the parameters option, *Selection?* greatly influence Divide behavior:
  - Image: Single: do all the steps for the image selected at step 1.
  - Image: Multi if CROP: do all the steps for the image selected at step 1, except if the image is a crop, in which case, do all the steps for all crop images detected in the folder containing the selected image.

For the following possibilities, select a folder instead of an image at step 1.

- Folder: CROPs only: do all the steps for all clean crops detected in the folder selected at step 1 and its subfolders.
- Folder: All but CROPs: do all the steps for all clean images (but not crops) detected in the folder selected at step 1 and its subfolders.
- Folder: All: do all the steps for all clean images detected in the folder selected at step 1 and its subfolders.

Image naming convention has to be followed for images to be found in the (sub)folders.

- If *Batch divide?* is ticked in the parameters, only step 1 is apparent to the users. Default settings will be used. It doesn't work with "Median" or "Mean \* coefficient" thresholds which requires user interaction.

## 4.4 Analyse Tool

1. Select Ratio image. If a RoiSet corresponding to the image is found, it is loaded in the ROI Manager.
2. Manage regions of interest then click "OK".




Figure 12 – Manage regions.

3. Choose to add another image or not. If yes, go back to step 1. If no, continue.






Figure 13 – Add stack?

4. Set parameters then click "OK". If *Plot Ratio vs intensity?* drop down menu is on *Vs intensity* or *Vs intensity ratio*, choose fluorescence image(s) and select region of interest for each Ratio image.




Figure 14 – Set parameters.

5. Four graphs are displayed: mean, standard deviation, normalized distribution of means and normalized distribution of standard deviation, versus time. If *Plot Ratio vs intensity?*

drop down menu was on *Vs intensity* or *Vs intensity ratio*, a dynamic 3D plot window is displayed as well in case of timelapses, or a 2D plot in case of single timepoints.




Figure 15 – Graphs.

## Notes

- It is possible to replot data after modifying regions of interest (Update or Delete in the ROI Manager), and/or deleting image, then clicking on Analyse tool.
- To add a region to an existing image, open it again with "Add another stack". It will select the correct window and allow you to add more regions.
- It is possible to load a previous analysis by selecting the text file generated by the analyse tool instead of a Ratio image (see File Management).

## 4.5 Param Tool

*Donor* and *Acceptor* fields are used to automatically detect acceptor when selecting donor, and to name images resulting from the different tools.

### Clean parameters:

Ticking *Batch clean?* will change the behavior of Clean tool. The user will choose a folder containing raw images instead of choosing an image. All pairs of Donor/Acceptor detected inside this folder will be processed without user interaction: automatic background removal and stack alignment if ticked.

*Align stacks?* defines default value during execution.

### Crop parameters:

Ticking *Batch crop?* will change the behavior of Crop tool. The user will choose a folder containing clean images instead of choosing an image. All pairs of Donor/Acceptor containing "\_clean" detected inside this folder and its subfolders will be processed.

If *Give a label to CROPs?* is ticked, a text file will be generated containing the crop numbers and associated labels.

### Divide parameters:

*Selection?* drop down menu will allow user to change Divide tool behavior:

- *Image: Single*: user selects clean image to be processed.

- *Image: Multi if CROP*: user selects clean image. If the image is a crop, all crops in the same folder will be processed. If it's not, the behavior is the same as Image: Single.

- *Folder: CROPs only*: user selects a folder. All clean crops contained in the folder and its subfolders will be processed.

- *Folder: All but CROPs*: user selects a folder. All clean images that are not crops contained in the folder and its subfolders will be processed.

- *Folder: All*: user selects a folder. All clean images contained in the folder and its subfolders will be processed. Image naming convention has to be followed for images to be found in the (sub)folders.

Ticking *Without confirmation?* will use *Threshold method*, *Coefficient*, *AutoLocalTh radius*, *Minimum threshold* and *Max range* values under without asking the user during execution of divide.

### Analyse parameters:

*Threshold method*, *Coefficient*, *AutoLocalTh radius*, *Minimum threshold*, *Min range* and *Max range* define default values during execution.

The threshold works as follow: for each image of donor stack, it computes the median of the area selected or mean \* coefficient of the area selected or global threshold or local threshold and put to 0 all pixels for which value is lower. After that, if *Minimum threshold* is higher than 0, set to 0 all pixels under that value.

Parameters

Donor: LUC

Acceptor: YFP

Clean parameters:

- ☐ Batch clean? (see Help)
- ☒ Align stacks?

Crop parameters:

- ☐ Batch crop? (see Help)
- ☒ Give a label to CROPs? (see Help)

Divide parameters:

Selection? (see Help): Image: Multi if CROP

- ☐ Without confirmation? (see Help)

Threshold method: AutoTh - Chastagnier

Coefficient: 1

AutoLocalTh radius: 10

Minimum threshold: 0

Min range: 0.200

Max range: 2

Analyse parameters:

Time between images (s): 30

Slice used to normalize: 5

- ☒ Save ROIs?

Plot Ratio vs intensity? No

Folder and data name: data

Help Cancel OK

*Time between images (s)*, *Slice used to normalize*, *Save ROIs?*, *Plot Ratio vs intensity?* and *Folder and data name* define default values during execution.

Setting *Plot Ratio vs intensity?* drop down menu on *Vs intensity* or *Vs intensity ratio* will allow the user to plot a 3D interactive window that displays the mean of the ROIs versus time and intensity (or intensity ratio) for timelapse or a 2D window for single timepoints.

## 4.6 Min Max Tool

1. Select the folder containing the "Ratio" images.
2. Select the range to apply to the set of images.

## 5 Output files structure

Following the naming convention, base images are named xxxDonor.tif and xxxAcceptor.tif. Where xxx is any character string, Donor and Acceptor are character strings from the parameters' fields Donor and Acceptor.

*Clean* process will create a subfolder xxx in which everything will be put, starting with clean images: xxxDonor\_clean.tif and xxxAcceptor\_clean.tif.

*Crop* process will extract a sub area of the clean images and save them as xxxCROP<i>\_Donor\_clean.tif and xxxCROP<i>\_Acceptor\_clean.tif, where <i> is crop's number. It will also save regions in cropAREAs.zip, and the crop's number with a label in file CROPs.txt if *Give a label to CROPs?* is ticked in parameters.

*Divide* process will create the ratiometric image xxxRatio.tif with threshold values used and/or method in text file thresholdUsed.txt. If CROP images are processed, the image and text file will be named xxxCROP<i>\_Ratio.tif and thresholdUsedCROP<i>.txt, where <i> is the crop's number.

*Analyse* process will use the character string entered in parameters' field *Folder and data name*, which will be referred as Data in the following. A text file Data.txt is created containing the ending part of the path to images used in the analysis. It allows the user to select it instead of an image during execution of *Analyse* to quickly load all the images contained in the file and their corresponding Regions of interest. A subfolder Data is created which will contain the following files: spreadsheets Data.csv and Data.xls containing mean, standard deviation and distributions of regions of interest, DataNorm.csv and DataNorm.xls containing the same data but normalized at a given slice, and Image<i>\_ROI<j>.png for each image <i> (or Image<i>\_ROI<j>-<k>.png if image contains multiple regions of interest).




Figure 16 – Files structure.

## 6 Changing default values

Each time the macro is loaded in Fiji, global variables are initialized to default values defined in the code. The user can then change them in the *Parameters* tool. If a value is always the same but differs from default, it might come in handy to modify it permanently instead of having to change it every time.

To do so, drag and drop the file `Luminescence_Ratio_Analyzer.ijm` onto Fiji. It will open the code editor. The highlighted section in Figure 17 is the part to edit. Lines have the form:

```
var varName = value; // comment
```

What you have to change is the value, between "=" and ";". There are three types of values:

- booleans for which value is either *true* or *false* (they appear in blue in the code editor).
- numbers (with a dot separator, in case it's not integers).
- strings: characters between double quotes (they appear in pink in the code editor).

```

1 // Luminescence Ratio Analyzer version 0.26 (2016-11-14)
2 // by Yan Chastagnier
3
4 // Put this macro into the folder Fiji.app/macros/toolsets/
5 // Select it with the >> symbol at the extrem right of the toolbar.
6 // Use the Action Tools to run the different parts of the process.
7
8 // Global parameters default values : change them to alter permanently default behaviour on toolset startup
9
10 var donorName = "LUC"; // Donor
11 var acceptorName = "YFP"; // Acceptor
12 var batchClean = false; // Batch clean?
13 var alignStacks = true; // Align stacks?
14 var batchCrop = false; // Batch crop?
15 var nameCROPs = true; // Give a label to CROPs?
16 var divideSelection = "Image: Multi if CROP";
17 var batchDivide = false; // Batch divide?
18 var thresholdMethod = "AutoTh - Chastagnier";
19 var coefMultMeanROI = 1; // Coefficient
20 var radiusLocalTh = 10; // AutoLocalTh radius
21 var overallMinThreshold = 0; // Minimum threshold
22 var rangeMin = 0.2; // Min range
23 var rangeMax = 2; // Max range
24 var time_between_images_sec = 30; // Time between images (s)
25 var normalizeSlice = 5; // Slice used to normalize
26 var saveROIs = true; // Save ROIs?
27 var plot3D = "No"; // Plot Ratio vs intensity?
28 var dataFolderName = "data"; // Folder and data name
29
30 // Tools definition
31 macro "Clean donor and acceptor images Action Tool - B03 C059 T060dC Ta60dl Td60de T3f0da Taf0dn" {cleanImages();}
32 macro "Crop clean donor and acceptor images Action Tool - B03 C059 T060dC Ta60dr T3f0do Tbf0dp" {cropImages();}
33 macro "Divide clean images to get Ratio image Action Tool - B03 C059 T050dD Ta50di Td50dv T2f0di T5f0dd Tcf0de" {makeRatioIm
34 macro "Ratio analysis Action Tool - B03 C059 T050dA T850dn Tf50da T0f0dl T3f0dy Taf0ds Tff0de" {makeRatioAnalysis();}
35 macro "Get and set image location and size Action Tool - B03 C059" {setLoc();}
36 macro "Close everything Action Tool - B03 C059 T060dC Ta60dl Td60do T3f0ds Taf0de" {closeImages();}

```

Figure 17 – Code to change default values.